



Roto-Dry

Forced Hot Air Laboratory Slide Dryer



Mopec The Better Product People

Call Us at: **1-800-362-8491** or Visit Us on the Web at www.mopec.com

EXTRA LARGE DRYING CHAMBER

The 13" diameter x 5" (33 cm x 12.7 cm) deep chamber is built to accommodate virtually any size horizontal or vertical slide staining rack. Digital keypad touch controls for easy operation. Large LED temperatures display to easily monitor delicate staining procedures.

BUILT IN ELECTRONIC TIMER

Roto-Dry comes equipped with a simple electronic programmable timer/alarm that will automatically sound when the slide drying process is complete. After the alarm sounds, the Roto-Dry will continue to operate to guarantee a constantly warm, and ready-to-use drying chamber.

ROTATING PLATFORM FORCED HOT-AIR

The unique rotating platform combined with temperature controlled forced hot air assures quick and even slide drying, eliminating hot or cold spots. Average drying time: 15-20 minutes @70°C.

LESS TECH TIME

The Roto-Dry allows for quick and easy placement of a slide rack through a large hinged opening. Close the lid, set the temperature and time and walk away. Manual rotation of racks is no longer necessary to assure quick, even drying.

DURABLE AND EASY TO CLEAN

Attractive stainless steel construction guarantees years of reliable service. All surfaces of the Roto-Dry will easily clean with any standard laboratory disinfectant.



"Dries more evenly and consistently than any other dryer I've experienced."

Histology Dept. Supervisor

Roto-Dry Specifications

Dimensions:	12-1/4" H x 15-3/4" W x 21-1/2" L (31.2 cm H X 40 CM w X 54.6 CM L)
Chamber Inside:	13" Round (33 cm)
Dimension:	5" Deep (12.7 cm)
Electrical:	115/1/60
Temperature Range:	Ambient to 80°C

Product No. BK200



For a complete catalog detailing our products and services, please call or visit our web site



21750 Coolidge Highway • Oak Park, MI 48237
Toll Free: 800-362-8491 • Tel: 248-291-2040 • Fax: 248-291-2050

www.mopec.com